

DuraScan *Alpha*

- ✓ Open System and Perfect Adaptation
- ✓ High Speed, High Accuracy
- ✓ Texture Scanning

erkdigital.com.tr



brochure

**Dimensions (mm)**

350 x 370 x 420

**Weight (kg)**

5

**Scanning Accuracy**

<10µm

**Light Source**

LED, Blue Light, ANSI Lumen

**Connection**

USB 3.0 B-Type

**Scan Range (mm)**

100 x 100 x 75

**Processing Method**

Phase-shifting Optical Triangulation

**Camera Resolution**

1.3 Megapixel x 2

High Speed, High Accuracy

With its hardware of the highest quality and optimized software, you can scan within seconds! In addition to its high scanning speed, DuraScan Alpha acquires perfect scans under 10 microns of accuracy, according to ISO 12836 Standard.



Open System and Perfect Adaptation

DuraScan Alpha offers completely system, and due to its ability of integration with exocad®, you can have faster and more efficient results in your works.



Texture Scanning

Texture scanning function offers perfect dentist-technician communication and ability to design ideal restorations.

